

## Shoulder Immobilizer – Universal



Figure 1



Figure 2



Figure 3



Figure 4



Figure 5



Figure 6

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

### Intended for single patient use.



**WARNING:** If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.

### APPLICATION INSTRUCTIONS

1. Unfold swathe and position contoured portion on front mid-section. Wrap swathe around waist and secure (Figure 1).
2. Thread narrow end of humeral strap through slit to form a loop. Slide humeral strap onto affected arm and position above elbow (Figure 2).
3. Pull evenly on both ends of humeral strap until properly fitted around arm. Secure wide end of strap on back of chest swathe and narrow end on front (Figure 3).
4. Wrap and secure cuff on wrist of affected arm. Secure cuff to chest swathe (Figure 4).
5. Secure the end of shoulder strap without a buckle to the back of the swathe (Figure 5).
6. Pull strap over unaffected shoulder and secure to the bottom of the wrist cuff. Use slide buckle to adjust length of shoulder strap (Figure 6).

**NOTE:** Please refer to sew-in label on device for care instructions.

Please visit our website at [www.birdcronin.com](http://www.birdcronin.com) for additional product and warranty information or contact a customer care representative at **800-328-1095**.

**NOTICE:** This device is only one element in the overall treatment program administered by a medical professional. While every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.

**NOT MADE WITH  
natural rubber latex**