

# Comfor™ Shoulder Abduction Pillow




Figure 1

## COMPONENTS

- A. Shoulder Abduction Pillow (45° and 70° angle)      C. Circumferential Strap  
B. Padded Shoulder Strap with Wrist Cuff              D. Humeral Cuff




Figure 2



Figure 3



Figure 4



Figure 5



Figure 6

**NOT MADE WITH  
natural rubber latex**

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

**This device is intended for single patient use.**



**WARNING:** If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.

## APPLICATION INSTRUCTIONS

1. Open package and separate components. (Figure 1)
2. Position the shoulder abduction pillow (A) under the arm at 45° or 70° as directed by healthcare provider. Elbow and forearm should rest comfortably on the pillow. (Figure 2)
3. Use the circumferential strap (C) to secure the pillow in place. Attach the strap to the pillow by wrapping around the patient's mid-section and back onto the pillow. (Figure 3)
4. Wrap humeral cuff (D) around the bicep and secure to the pillow. (Figure 4)
5. Wrap wrist cuff, that is attached to the shoulder strap (B), around wrist of the affected arm. Attach hook on wrist cuff to the pillow. (Figure 5)
6. Bring padded shoulder strap (B) up and over unaffected shoulder and secure to the back of abduction pillow. (Figure 6)
7. Secure shoulder strap on the front by threading the strap through the buckle on wrist cuff. (Figure 6)
8. Adjust length of shoulder strap, placement of foam neck pad, and circumferential strap around abduction pillow if necessary.

**NOTE:** Please refer to sew-in label on device for care instructions.

Please visit our website at [www.birdcronin.com](http://www.birdcronin.com) for additional product information or contact a customer care representative at **800-328-1095**.

**NOTICE:** Every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort with this product, however, there is no guarantee that injury will be prevented with its use.