

B-Cool® Teri Shoulder Wrap and Kits

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

Intended for single patient use.



WARNING: This product should be used only under the supervision of a medical professional. If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.



CAUTION: Never apply a frozen or heated gel pack directly onto the skin.



CAUTION: Supervision is strongly advised when applying this product on children and geriatric patients. Patients with circulatory problems should consult a physician or healthcare professional prior to any treatment.

APPLICATION INSTRUCTIONS

NOTE: Before use, follow directions on B-Cool® Gel Packs for Cold or Heat therapy.

1. Open box, remove wrap and gel packs (if included with the kit).
2. Heat or freeze gel packs according to physicians instructions following directions on B-Cool® Gel Pack packaging.
3. Insert prepared gel packs into wrap pockets.
2. Release buckle and place wrap over affected shoulder - Teri cloth side down.
3. Align the center area of the wrap over the joint of the affected shoulder with the longer sections going across the upper shoulder/neck.
3. Position wrap so one pack lies over front of the chest and the other lies over the back.
4. Wrap strap around back and secure with front buckle.
3. Adjust straps as needed for comfort.

NOTE: Please refer to liner sew-in label for care instructions.

Please visit our website at www.birdcronin.com for additional product and warranty information or contact a customer care representative at **800-328-1095**.

NOTICE: This device is only one element in the overall treatment program administered by a medical professional. While every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.

